

UK Office

Everest Biotech Ltd

Cherwell Innovation Centre 77 Heyford Park Upper Heyford Oxfordshire OX25 5HD

Enquiries:

info@everestbiotech.com

Sales:

UK

sales@everestbiotech.com

Tech support:

support@everestbiotech.com

Tel: +44 (0)1869 238326

www.everestbiotech.com

Research Use Only. Not for diagnostic or therapeutic use.

EB06327 - Goat Anti-PITPNM / PITPNM1 Antibody

Size: 100µg specific antibody in 200µl



Target Protein

Principal Names: PITPNM, NIR2, DRES9, phosphatidylinositol transfer protein, membrane-associated, Drosophila retinal degeneration B, PYK2 N-terminal domain-interacting receptor 2, PITPNM1, FLJ44997, RDGB, RDGB1, RDGBA, RDGBA1, Rd9, phosphatidylinositol transfer protein, membrane-associated, retinal degeneration B

Official Symbol: PITPNM1

Accession Number(s): NP_004901.2; NP_001124320.1

Human GeneID(s): 9600

Non-Human GeneID(s): 18739 (mouse)

Important Comments: This antibody is expected to recognise both reported isoforms.

Immunogen

alpha 1

Peptide with sequence C-KARSISLKLDSEE, from the C Terminus of the protein sequence according to NP_004901.2; NP_001124320.1.

Please note the peptide is available for sale.

Purification and Storage

Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin

Aliquot and store at -20°C. Minimize freezing and thawing.

Applications Tested

Peptide ELISA: antibody detection limit dilution 1:8000.

Immunofluorescence: Strong expression of the protein seen in the cytoplasm and endoplasmic reticulum of A549 cells. Recommended concentration: 10µg/ml.

Flow Cytometry: Flow cytometric analysis of HepG2 cells. Recommended concentration: 10ug/ml.

Species Reactivity

Tested: Human

Expected from sequence similarity: Human, Mouse, Dog, Cow

EB06327 Immunofluorescence analysis of paraformaldehyde fixed A549 cells, permeabilized with 0.15% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml), showing cytoplasmic and Endoplasmic Reticulum staining. The nuclear stain is DAPI (blue). Negative control: Unimmunized goat IgG (10ug/ml) followed by Alexa Fluor 488 secondary antibody (2ug/ml).

EB06327 Flow cytometric analysis of paraformaldehyde fixed A549 cells (blue line), permeabilized with 0.5% Triton. Primary incubation 1hr (10ug/ml) followed by Alexa Fluor 488 secondary antibody (1ug/ml). IgG control:

Unimmunized goat IgG (black line) followed by Alexa Fluor 488 secondary antibody.